



More information on the website  
[mirror.radwag.com/de/info,w1,6QU](https://mirror.radwag.com/de/info,w1,6QU)

# Prüfgewicht der M1 Klasse - Plättchenform - 500 mg

OP-403-0075



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Konstruktion	
Material	Neusilber
Nominalwert	500 mg

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself