



More information on the website  
[mirror.radwag.com/en/info,w1,6YQ](https://mirror.radwag.com/en/info,w1,6YQ)

# Knob weights, Set (1 g - 1 kg), Wooden Box

OK-101-0034



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | E2              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 2 110 g         |
| Nominal value           | 1 g – 1 kg      |
| Quantity                | 13 pcs          |

## Compatible with (Additional Fee)

Tweezers

Brushes