



More information on the website  
[mirror.radwag.com/en/info,w1,6YQ](https://mirror.radwag.com/en/info,w1,6YQ)

# Knob weights, Set (1 g - 1 kg), Wooden Box

OK-101-0034



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	2 110 g
Nominal value	1 g – 1 kg
Quantity	13 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes