



More information on the website  
[mirror.radwag.com/de/info,w1,6Z5](http://mirror.radwag.com/de/info,w1,6Z5)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 100 g), Kunststoff-Etui

OK-402-0072



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	211,11 g
Nominalwert	1 mg – 100 g
Anzahl	21 St.

## Kompatibel mit (Additional Fee)

Pinsel