



More information on the website  
[mirror.radwag.com/en/info,w1,6Z5](https://mirror.radwag.com/en/info,w1,6Z5)

# M1 Mass Standard, Knob Weights, Set (1 mg - 100 g), Plastic Box

OK-402-0072



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	M1
Weight	211.11 g
Nominal value	1 mg – 100 g
Quantity	21 pcs

## Compatible with (Additional Fee)

Brushes