



More information on the website  
[mirror.radwag.com/en/info,w1,74S](https://mirror.radwag.com/en/info,w1,74S)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 10 kg), Plastic Box

OK-202-0061



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 21 110 g        |
| Nominal value           | 1 g – 10 kg     |
| Quantity                | 17 pcs          |

## Compatible with (Additional Fee)

Tweezers

Brushes