



More information on the website
mirror.radwag.com/de/info,w1,7TT

Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 500 g), Holz- Etui

OK-401-0082



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	1 111,11 g
Nominalwert	1 mg - 500 g
Anzahl	24 St.

Kompatibel mit (Additional Fee)

Pinsel