



More information on the website
mirror.radwag.com/tr/info,w1,8F2

KIT 85 Yoğunluk Tespit Kiti



The drawings, photos and graphics used are for illustrative purposes only.

Genel Açıklama

KIT 85 is a density determination kit intended for laboratory balances featuring a weighing pan with $\varnothing 85$ and $\varnothing 100$ mm.

This stainless steel accessory allows determining density of:

solids (including granulates); liquids.

KIT 85 features an additional set of weighing pans for testing granulate of density higher, and lower than 1 g/cm³.

İle uyumlu

AS X7 Analitik Teraziler
AS X2 PLUS Analitik Teraziler
XA 5Y.A Analitik Terazi
AS 5Y Analitik Teraziler

Filtreleri için Analitik Terazi XA 5Y.F
XA 5Y Analitik Terazi
AS R2 PLUS Analitik Teraziler
AS R1 PLUS Analitik Teraziler


Cihaz boyutları




KIT 85




KIT 128



KIT 195



KIT 195.M



Granule density
determination KIT