



More information on the website  
[mirror.radwag.com/it/info,w1,8MB](https://mirror.radwag.com/it/info,w1,8MB)

# Cilindrico, set (1 g - 5 kg), imballaggi in plastica

OK-302-0030



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Classe OIML	F2
Costruzione	
Materiale	Stainless steel
Peso	11 110 g
Valore nominale	1 g – 5 kg
Numero	16 pcs

## Compatibile con (Additional Fee)

Pinzette

Pennelli