



More information on the website
mirror.radwag.com/en/info,w1,8V1

F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 500 g), Plastic Box

OK-202-0057



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	1 110 g
Nominal value	1 – 500 g
Quantity	12 pcs

Compatible with (Additional Fee)

Tweezers

Brushes