



More information on the website
mirror.radwag.com/en/info,w1,93J

F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 5 kg), Wooden Box

OK-201-0049



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	11 110 g
Nominal value	1 g - 5 kg
Quantity	16 pcs

Compatible with (Additional Fee)

Tweezers

Brushes