



More information on the website
mirror.radwag.com/de/info,w1,98E

Prüfgewicht der E2 Klasse - Plättchenform - 200 mg

OP-103-0069



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

| Messtechnische Parameter | |
|--------------------------|-----------|
| OIML-Klasse | E2 |
| Konstruktion | |
| Material | Neusilber |
| Nominalwert | 200 mg |

Kompatibel mit (Additional Fee)

Pinzetten

Pinself