



More information on the website  
[mirror.radwag.com/en/info,w1,9CJ](https://mirror.radwag.com/en/info,w1,9CJ)

## BT 1000 µl Tips - 250 pcs

PT-112-0001



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

OMNITIP™ BT 1000 µl tips for pipettes are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

**Capacity: 1000 µl Colour: blue Set contains 250 tips in a plastic bag with a zipper closure**

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

## Compatible with (Additional Fee)

