



More information on the website  
[mirror.radwag.com/en/info,w1,9CJ](http://mirror.radwag.com/en/info,w1,9CJ)

## BT 1000 µl Tips - 250 pcs

PT-112-0001



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

OMNITIP™ BT 1000 µl tips for pipettes are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

Capacity: 1000 µl Colour: blue Set contains 250 tips in a plastic bag with a zipper closure

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

## Compatible with (Additional Fee)

Workstation for Pipettes Calibration