



More information on the website  
[mirror.radwag.com/en/info,w1,9GC](https://mirror.radwag.com/en/info,w1,9GC)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 g - 10 kg), Wooden Box

OK-201-0068



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	21 110 g
Nominal value	1 g – 10 kg
Quantity	17 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes