



More information on the website  
[mirror.radwag.com/de/info,w1,9KC](http://mirror.radwag.com/de/info,w1,9KC)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 200 g), Holz- Etui

OK-401-0081



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	611,11 g
Nominalwert	1 mg – 200 g
Anzahl	23 St.

## Kompatibel mit (Additional Fee)

Pinsel