



More information on the website
mirror.radwag.com/fr/info,w1,9VJ

Boîte anti-courant d'air pour microbalances

WX-001-0117



The drawings, photos and graphics used are for illustrative purposes only.





Paramètres Techniques

Paramètres physiques	
Dimension	560x350x252 mm
Dimensions de colis L x P x H	695x485x385 mm
Masse nette	7,6 kg
Masse brute	10,9 kg


Compatible avec (Additional Fee)

Microbalance pour étalonnage des pipettes XA 5Y.M.A.P
 Microbalance pour étalonnage des pipettes MYA 5Y.P
 Microbalance pour filtres MYA 5Y.F
 Microbalance XA 5Y.M.A


Microbalance MYA 5Y
 Microbalance XA 5Y.M
 Microbalance pour stents XA 5Y.M.A.S






Dimensions d'appareil L x P x H




[Anti-draft chamber for balances with 128x128 mm weighing pan](#)




[Anti-draft chamber for X2 and R2 precision balances](#)




[Anti-draft chamber for microbalances](#)




[Draft shield for precision balances](#)




[Anti-draft chamber for APP-00-SK mass comparators](#)



[Anti-draft chamber for APP-KO-SK mass comparators](#)



[Anti-draft chamber for WAY-KO mass comparators](#)



[Anti-draft chamber for XA 5Y / XA 5Y.A balances](#)