



More information on the website  
[mirror.radwag.com/pl/info,w1,9VJ](https://mirror.radwag.com/pl/info,w1,9VJ)

## Szafka przeciwpodmuchowa do mikrowag

WX-001-0117



Użyte rysunki, zdjęcia, grafiki mają charakter poglądowy.


### Dane techniczne

| Parametry fizyczne           |                |
|------------------------------|----------------|
| Wymiar                       | 560x350x252 mm |
| Wymiary opakowania S x G x W | 695x485x385 mm |
| Masa netto                   | 7,6 kg         |
| Masa brutto                  | 10,9 kg        |


### Pasuje do (Dodatkowo płatne)

Mikrowaga do kalibracji pipet XA 5Y.M.A.P  
Mikrowaga do kalibracji pipet MYA 5Y.P  
Mikrowaga do filtrów MYA 5Y.F  
Mikrowaga XA 5Y.M.A

Mikrowaga MYA 5Y  
Mikrowaga XA 5Y.M  
Mikrowaga do ważenia stentów XA 5Y.M.A.S




## Wymiary urządzenia S x G x W




[Anti-draft chamber for balances with 128x128 mm weighing pan](#)


[Anti-draft chamber for X2 and R2 precision balances](#)




[Anti-draft chamber for microbalances](#)




[Draft shield for precision balances](#)




[Anti-draft chamber for APP-00-SK mass comparators](#)



[Anti-draft chamber for APP-KO-SK mass comparators](#)



[Anti-draft chamber for WAY-KO mass comparators](#)



[Anti-draft chamber for XA 5Y / XA 5Y.A balances](#)