



More information on the website
mirror.radwag.com/us/info,w1,9VJ

Anti-Draft Chamber for Microbalances

WX-001-0117



The drawings, photos and graphics used are for illustrative purposes only.


Datasheet

Physical parameters	
Dimension	560x350x252 mm
Packaging dimensions W x D x H	695x485x385 mm
Net weight	7,6 kg
Gross weight	10,9 kg


Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 MYA 5Y.F Microbalance for Filters
 XA 5Y.M.A Microbalance

MYA 5Y Microbalance
 XA 5Y.M Microbalance
 XA 5Y.M.A.S Microbalance for stents weighing




Device dimensions W x D x H




[Anti-draft chamber for balances with 128x128 mm weighing pan](#)


[Anti-draft chamber for X2 and R2 precision balances](#)




[Anti-draft chamber for microbalances](#)




[Draft shield for precision balances](#)




[Anti-draft chamber for APP-00-SK mass comparators](#)



[Anti-draft chamber for APP-KO-SK mass comparators](#)



[Anti-draft chamber for WAY-KO mass comparators](#)



[Anti-draft chamber for XA 5Y / XA 5Y.A balances](#)