



More information on the website
mirror.radwag.com/it/info,w1,9XE

Pipetta automatica 100 μ l ÷ 1000 μ l - RW8-106-20-9

PA-203-0006



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri metrologici

Volume	100 – 1000 μ l
Systematic error	2,0% / 0,6%
Random error	0,7% / 0,2%
Tips	1000 μ l



Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P
 Microbilancia per la calibrazione delle pipette MYA 5Y.P

Bilance automatiche per la calibrazione delle pipette