



radwag.com



## 100 µl ÷ 1000 µl - RW8-106-20-9 Otomatik Pipetler

More information on the website  
[mirror.radwag.com/tr/info/w1,9XE](http://mirror.radwag.com/tr/info/w1,9XE)



The drawings, photos and graphics used are for illustrative purposes only.

### Teknik Veriler

#### Metrolojik parametreler

Hacim	100 µl ÷ 1000 µl
Systematic error	2,0% / 0,6%
Tekrarlanabilirlik error	0,7% / 0,2%
Uçlar	1000 µl



### İle uyumlu

XA 5Y.M.A.P Mikro Teraziler  
MYA 5Y.P pipetlerini kalibre etmek için Mikro Terazi

Pipet kalibrasyonu için otomatik teraziler