



More information on the website
mirror.radwag.com/en/info,w1,9XL

F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 mg - 1 kg), Plastic Box

OK-202-0051



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	2 111.11 g
Nominal value	1 mg – 1 kg
Quantity	25 pcs

Compatible with (Additional Fee)

Tweezers

Brushes