



# Antivibration Table for pipette calibration

WX-001-0481

More information on the website  
[mirror.radwag.com/en/info,w1,A4X](http://mirror.radwag.com/en/info,w1,A4X)



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Construction	mild steel, particle board, granite
Load-bearing capacity	50 kg
Physical parameters	
Packaging dimensions W x D x H	1.6×0.9×1 m
Net weight	100 kg
Gross weight	150 kg
Length	1500 mm
Width	800 mm
Height	800±100 mm

## Compatible with (Additional Fee)

Automatic Balances for Pipettes Calibration

# Device dimensions W x D x H

