



Antivibration Table for pipette calibration

WX-001-0481

More information on the website
mirror.radwag.com/us/info,w1,A4X



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	mild steel, particle board, granite
Load-bearing capacity	50 kg
Physical parameters	
Packaging dimensions W x D x H	1,6×0,9×1 m
Net weight	100 kg
Gross weight	150 kg
Length	1500 mm
Width	800 mm
Height	800±100 mm



Compatible with (Additional Fee)

Device dimensions W x D x H

