



More information on the website
mirror.radwag.com/us/info,w1,A4X

Laboratory Antivibration Table for pipette calibration



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	61 kg
Gross weight	65-70 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with (Additional Fee)

Automatic Balances for. Pipettes Calibration