



More information on the website  
[mirror.radwag.com/it/info,w1,A7T](https://mirror.radwag.com/it/info,w1,A7T)

# Pipetta automatica 2 µl ÷ 20 µl - RW8-102-20-9

PA-203-0002



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

Volume	2 – 20 µl
Systematic error	3,0% / 0,9%
Random error	2,0% / 0,4%
Tips	300 µl



## Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P  
 Microbilancia per la calibrazione delle pipette MYA 5Y.P

Bilance automatiche per la calibrazione delle pipette  
 Puntali per pipette