



More information on the website  
[mirror.radwag.com/en/info,w1,AGB](https://mirror.radwag.com/en/info,w1,AGB)

# Knob weights, Set (1 g - 50 g), Wooden Box

OK-101-0030



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	110 g
Nominal value	1 – 50 g
Quantity	8 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes