



More information on the website
mirror.radwag.com/us/info,w1,A07

PIPETTES Extended Module (5Y)

WX-016-0331



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PIPETTES Extended Module is designed for checking and calibrating single-channel pipettes according to ISO 8655. We have implemented the module in the 5Y terminal as an optional feature that you can activate. If you choose to do so, you will receive a lifetime license. The PIPETTES Extended Module creates its own web page (works together with the Web Editor application) and is based on the PIPETTES software. The Web Editor is used to manage databases, edit pipette calibration orders and generate reports on calibration with the use of a web browser.

Compatible with (Additional Fee)

XA 5Y.M.A Microbalance
XA 5Y Analytical Balance
XA 5Y.A Analytical Balance

XA 5Y.M Microbalance
MYA 5Y Microbalance