



# Knob weights, Set (1 g - 50 g), Plastic Box

OK-502-0004

More information on the website  
[mirror.radwag.com/en/info,w1,B0A](http://mirror.radwag.com/en/info,w1,B0A)



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E1
Construction	
Material	Stainless steel
Weight	110 g
Nominal value	1 – 50 g
Quantity	8 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes