



More information on the website
mirror.radwag.com/en/info,w1,B4I

F2 Mass Standard, Flat polygonal sheets, Set (1 mg - 500 mg), Wooden Box

OK-301-0026



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|---------------|
| OIML Class | F2 |
| Construction | |
| Material | Nickel silver |
| Weight | 1.11 g |
| Nominal value | 1 mg - 500 mg |
| Weights quantity | 12 |

Compatible with (Additional Fee)

Tweezers

Brushes