



More information on the website
mirror.radwag.com/it/info,w1,B5C

Tavolo di pesata professionale PSW

WX-001-0020



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri fisici	
Componenti del kit	table for calibration, anti-vibration table
Dimensioni del pacco L x P x A	1160×820×986 mm
Peso net	150 kg (table 110 + stone 40)
Peso lordo	171 kg
Lunghezza	1100 mm
Larghezza	750 mm
Altezza	1100 mm
Number of workstations	1



Compatibile con (Additional Fee)





Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P
 Microbilancia per la calibrazione delle pipette MYA 5Y.P
 Bilancia di precisione 5Y.PM

Bilancia di precisione CY10
 Bilancia analitica per filtri XA 5Y.F
 Microbilancia MYA 5Y

Microbilancia per filtri MYA 5Y.F
 AS X7 Bilancia analitica
 Microbilancia XA 5Y.M.A
 Bilancia di precisione C32.PM
 Bilancia analitica XA 5Y.A
 UYA 5Y Ultra-Microbilancia
 Bilancia analitica AS 5Y
 Bilancia di precisione PS 5Y
 Bilance PS 5Y per controllo CPC

Bilancia analitica XA 5Y
 Bilancia analitica AS R2 PLUS
 Microbilancia XA 5Y.M
 UYA 5Y.KO Comparatore di massa manuale
 Bilancia analitica AS R1 PLUS
 UYA 5Y.F Ultra-microbilancia per filtri
 Microbilancia per stent XA 5Y.M.A.S
 Moduli MAS.1.Y

Device dimensions L x P x A



—SDKP

—SDKP DUAL