



More information on the website  
[mirror.radwag.com/tr/info,w1,B97](https://mirror.radwag.com/tr/info,w1,B97)

## KIT 128 Yoğunluk Tespit Kiti



The drawings, photos and graphics used are for illustrative purposes only.

### Genel Açıklama

KIT 128 is a density determination kit intended for laboratory balances featuring 128x128 weighing pan.

This stainless steel accessory allows determining density of:

solids, liquids.


For testing granulates of density value higher or lower than 1 g/cm<sup>3</sup> use KIT 85.

### İle uyumlu

PUK için PS 5Y terazileri  
PS X7 Hassas Teraziler  
PS 5Y Hassas Teraziler  
PS X2 Hassas Teraziler

PS R1 Hassas Teraziler  
PUK için PS X7 terazileri  
PS R2 Hassas Teraziler  
PUK için PS X2 terazileri


# Cihaz boyutları




KIT 85




KIT 128



KIT 195



KIT 195.M



Granulte density  
determination KIT