



More information on the website
mirror.radwag.com/us/info,w1,BK3

DJ-05 Antistatic ionizer

WX-001-0098



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet


Power supply	100 – 240 V AC 50/60 Hz
Operating temperature	0 – +40 °C
Relative humidity	5% – 80%
Distance form a sample	Ok. 5 – 50 mm
Ion source	stainless steel
Physical parameters	
Frame dimensions	175×70×172 mm
Packaging dimensions W x D x H	375×375×130 mm
Net weight	1,3 kg
Gross weight	2,5 kg

Compatible with (Additional Fee)


MYA 5Y.P Microbalance for pipettes calibration
 MYA 5Y.F Microbalance for Filters

MYA 5Y Microbalance
 UYA 5Y.F Ultra-Microbalance for Filters


Device dimensions W x D x H



DJ-04



DJ-05



power supply
