



More information on the website  
[mirror.radwag.com/en/info,w1,BQI](https://mirror.radwag.com/en/info,w1,BQI)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 kg - 5 kg), Wooden Box

OK-301-0035



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	10 000 g
Nominal value	1 – 5 kg
Quantity	4 pcs

## Compatible with (Additional Fee)

Brushes