



More information on the website  
[mirror.radwag.com/us/info,w1,BS0](http://mirror.radwag.com/us/info,w1,BS0)

# Hybrid Protective Glass 5Y, CY10

WX-001-0435



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

Protective cover is made of transparent plastic, secures a display against scratches, dust etc. It is intended for balances X2 and R2 type.

Set - 10 pcs.

## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration  
 UMA 5Y.F Automatic Weighing System  
 5Y.PM Precision balance  
 MYA 5Y.F Microbalance for Filters  
 XA 5Y.M.A Microbalance  
 CY10 Scales for PGC  
 PUE CY10 Weighing Terminal  
 CY10 Multifunctional Scale  
 Automatic Balances for. Pipettes Calibration  
 XA 5Y.A Analytical Balance  
 WAY 5Y.KO Manual Mass Comparator  
 UYA 5Y Ultra-Microbalance  
 AS 5Y Analytical Balance  
 AK-4 Automatic Mass Comparator  
 PS 5Y Precision Balance

APP 5Y.KO Manual Mass Comparator  
 Robotic Mass Comparators  
 Density measurement  
 MYA 5Y Microbalance  
 XA 5Y Analytical Balance  
 Manual susceptometers - Magnetization  
 5Y.PM.KB Manual Mass Comparator  
 HRP 5Y.KO Mass Comparator  
 Automatic Vacuum Mass Comparators  
 PMV 5Y Moisture Analyzer  
 Automatic susceptometers - Magnetization  
 XA 5Y.KO Manual Mass Comparator  
 HRP 5Y.KB Mass Comparator  
 XA 5Y.M Microbalance  
 RMC F Robotic Weighing System  
 RB F Robotic Weighing System

PS 5Y Balances for PGC  
AKM-2 Automatic Mass Comparator  
CY10 Precision balance  
Analytical Balance for Filters XA 5Y.F  
UMA Automatic Mass Comparator

UYA 5Y.KO Manual Mass Comparator  
UYA 5Y.F Ultra-Microbalance for Filters  
Nano-Comparator NANO.AK-4.500.5Y  
XA 5Y.M.A.S Microbalance for stents weighing  
AK-6 5Y.F Automatic Weighing System