



More information on the website  
[mirror.radwag.com/en/info,w1,BSB](http://mirror.radwag.com/en/info,w1,BSB)

# LS2208 Barcode Scanner

WX-006-0005



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

LS2208 barcode scanner is a manual device designed for operation in adverse industrial conditions. It integrates with RADWAG terminals, PUE H315, PUE 5, PUE 7.1, HY10, PUE C32, and balances, 5Y, R2, X2, H7.

LS2208 is a high quality handheld barcode scanner. It has been designed for intense and continuous operation in retail sales, libraries, health centers and pharmacy, etc. Its light and ergonomic casing provides comfort of use. The function of dynamic scanning and wide scanning range, enable LS2208 to collect barcodes in a quick and reliable way, increasing work pace and efficiency. Intellistand™ system provides scanning of codes placed in front of the scanner, and easy switching between handheld and handfree mode (operation on a stand). A patented design of a single printed-circuit board eliminates failure endangered parts, i.e. connections between multiple boards. Scanning part of LS2208 with liquid polymere protects its module from shock based damages. LS2208 scanner features like integrated operation of various systems, universal cable or new symbols detection guarantee indispensable protection of an investment now and in the future.

## Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer  
 HX7 H Stainless steel Multifunctional Scale

HX7.4P2 Multifunctional Beam Scales  
 WAY 5Y.KO Manual Mass Comparator

HX7.4N H Multifunctional Stainless Steel Ramp Scale	UYA 5Y Ultra-Microbalance
HX7.4C Multifunctional Platform Scale	MA 5Y.IC Moisture Analyzer
HX7 HR Multifunctional Scale	MA R Moisture Analyzer
XA 5Y.M.A.P Microbalance for pipettes calibration	AS 5Y Analytical Balance
WPY Multifunctional Scale	AK-4 Automatic Mass Comparator
MYA 5Y.P Microbalance for pipettes calibration	PS 5Y Balances for PGC
C32 Multifunctional Scale	AKM-2 Automatic Mass Comparator
HY10.HRP.H Stainless steel High Resolution Scale	PS R2.H Precision Balance
MA 5Y Moisture Analyzer	HY10.H Scales for PGC
5Y.PM Precision Balance	CY10 Precision Balance
MYA 5Y.F Microbalance for Filters	XA 5Y.F Analytical Balance for Filters
HX7.4P2 H Multifunctional Stainless Steel Beam Scales	WLY Precision Balance
TMX H Stainless steel Multifunctional Scale	UMA Automatic Mass Comparator
WLY Scales for PGC	APP 5Y.KO Manual Mass Comparator
AS X7 Analytical Balance	Robotic Mass Comparators
PS X7 Balances for PGC	Density Measurement
TMX Multifunctional Scale	MYA 5Y Microbalance
XA 5Y.M.A Microbalance	PS R1 Precision Balance
PS R2 Precision Balance	XA 5Y Analytical Balance
WLC X7 Precision Balance	PUE HX7 Weighing Terminal
HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale	5Y.PM.KB Manual Mass Comparator
CY10 Scales for PGC	AS R2 PLUS Analytical Balance
DWT/RC/HYF Checkweighers	PUE C32 Weighing Terminal
PUE CY10 Weighing Terminal	PUE 7.1 Weighing Terminal
CY10 Multifunctional Scale	PUE 5 Weighing Terminal
HX7.4P H Multifunctional Stainless Steel Pallet Scales	HRP 5Y.KO Mass Comparator
HY10 Multifunctional Scale	PS 5Y Precision Balance
HY10 H Stainless Steel Multifunctional Scale	Automatic Vacuum Mass Comparators
HY10.HRP High Resolution Scale	PMV 5Y Moisture Analyzer
C32.PM Precision Balance	XA 5Y.KO Manual Mass Comparator
PUE 7.1.HRP High Resolution Scale	HRP 5Y.KB Mass Comparator
Automatic Balances for Pipettes Calibration	XA 5Y.M Microbalance
HX7 Multifunctional Scale	UYA 5Y.KO Manual Mass Comparator
MA X7.A Moisture Analyzer	AS R1 PLUS Analytical Balance
HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version	UYA 5Y.F Ultra-Microbalance for Filters
PS X7 Precision Balance	H315 H Waterproof Platform Scale
WPY Scales for PGC	HX7.4P Multifunctional pallet scales
XA 5Y.A Analytical Balance	Nano-Comparator NANO.AK-4.500.5Y
HX7.4 H Multifunctional Stainless Steel Platform Scale	XA 5Y.M.A.S Microbalance for stents weighing