



More information on the website
mirror.radwag.com/en/info,w1,BSB

LS2208 Barcode Scanner

WX-006-0005



The drawings, photos and graphics used are for illustrative purposes only.

Overview

LS2208 barcode scanner is a manual device designed for operation in adverse industrial conditions. It integrates with RADWAG terminals, PUE H315, PUE 5, PUE 7.1, HY10, PUE C32, and balances, 5Y, R2, X2, H7.

LS2208 is a high quality handheld barcode scanner. It has been designed for intense and continuous operation in retail sales, libraries, health centers and pharmacy, etc. Its light and ergonomic casing provides comfort of use. The function of dynamic scanning and wide scanning range, enable LS2208 to collect barcodes in a quick and reliable way, increasing work pace and efficiency. Intellistand™ system provides scanning of codes placed in front of the scanner, and easy switching between handheld and handfree mode (operation on a stand). A patented design of a single printed-circuit board eliminates failure endangered parts, i.e. connections between multiple boards. Scanning part of LS2208 with liquid polymere protects its module from shock based damages. LS2208 scanner features like integrated operation of various systems, universal cable or new symbols detection guarantee indispensable protection of an investment now and in the future.

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
MUYA Modules

WPY Scales for PGC
XA 5Y.A Analytical Balance

HX7 H Stainless steel Multifunctional Scale
HX7.4N H Multifunctional Stainless Steel Ramp Scale
HX7.4C Multifunctional Platform Scale
HX7 HR Multifunctional Scale
XA 5Y.M.A.P Microbalance for pipettes calibration
WPY Multifunctional Scale
MYA 5Y.P Microbalance for pipettes calibration
C32 Multifunctional Scale
HY10.HRP.H Stainless steel High Resolution Scale
MA 5Y Moisture Analyzer
5Y.PM Precision Balance
MYA 5Y.F Microbalance for Filters
HX7.4P2 H Multifunctional Stainless Steel Beam Scales
TMX H Stainless steel Multifunctional Scale
WLY Scales for PGC
MA X2.IC.A Moisture Analyzers
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balances
PS X7 Balances for PGC
TMX Multifunctional Scale
XA 5Y.M.A Microbalance
PS R2 Precision Balance
WLC X7 Precision Balance
HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale
CY10 Scales for PGC
PS X2 Balances for PGC
DWT/RC/HYF Checkweighers
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
HX7.4P H Multifunctional Stainless Steel Pallet Scales
MA X2.A Moisture Analyzers
HY10 Multifunctional Scale
HY10 H Stainless Steel Multifunctional Scale
HY10.HRP High Resolution Scale
C32.PM Precision Balance
PUE 7.1.HRP High Resolution Scale
Automatic Balances for. Pipettes Calibration
HX7 Multifunctional Scale
MA X7.A Moisture Analyzer
HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version
PS X7 Precision Balance

HX7.4 H Multifunctional Stainless Steel Platform Scale
HX7.4P2 Multifunctional Beam Scales
WAY 5Y.KO Manual Mass Comparator
UYA 5Y Ultra-Microbalance
MA 5Y.IC Moisture Analyzer
MA R Moisture Analyzer
AS 5Y Analytical Balance
AK-4 Automatic Mass Comparator
PS 5Y Balances for PGC
AKM-2 Automatic Mass Comparator
PS R2.H Precision Balance
HY10.H Scales for PGC
CY10 Precision Balance
XA 5Y.F Analytical Balance for Filters
WLY Precision Balance
UMA Automatic Mass Comparator
APP 5Y.KO Manual Mass Comparator
Robotic Mass Comparators
Density Measurement
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
5Y.PM.KB Manual Mass Comparator
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PUE HX7 Weighing Terminal
PUE 5 Weighing Terminal
HRP 5Y.KO Mass Comparator
PS 5Y Precision Balance
Automatic Vacuum Mass Comparators
PMV 5Y Moisture Analyzer
XA 5Y.KO Manual Mass Comparator
HRP 5Y.KB Mass Comparator
XA 5Y.M Microbalance
UYA 5Y.KO Manual Mass Comparator
AS R1 PLUS Analytical Balance
UYA 5Y.F Ultra-Microbalance for Filters
H315 H Waterproof Platform Scale
HX7.4P Multifunctional pallet scales
Nano-Comparator NANO.AK-4.500.5Y
XA 5Y.M.A.S Microbalance for stents weighing