



More information on the website  
[mirror.radwag.com/us/info,w1,BSB](https://mirror.radwag.com/us/info,w1,BSB)

# LS2208 Barcode Scanner

WX-006-0005



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

LS2208 barcode scanner is a manual device designed for operation in adverse industrial conditions. It integrates with RADWAG terminals, PUE H315, PUE 5, PUE 7.1, HY10, PUE C32, and balances, 5Y, R2, X2, H7.

LS2208 is a high quality handheld barcode scanner. It has been designed for intense and continuous operation in retail sales, libraries, health centers and pharmacy, etc. Its light and ergonomic casing provides comfort of use. The function of dynamic scanning and wide scanning range, enable LS2208 to collect barcodes in a quick and reliable way, increasing work pace and efficiency. Intellistand™ system provides scanning of codes placed in front of the scanner, and easy switching between handheld and handfree mode (operation on a stand). A patented design of a single printed-circuit board eliminates failure endangered parts, i.e. connections between multiple boards. Scanning part of LS2208 with liquid polymere protects its module from shock based damages. LS2208 scanner features like integrated operation of various systems, universal cable or new symbols detection guarantee indispensable protection of an investment now and in the future.

## Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer  
HX7 H Multifunctional Scale

HX7.4P2 Multifunctional Beam Scales  
WAY 5Y.KO Manual Mass Comparator

HX7.4N H Multifunctional Stainless Steel Ramp Scale  
 HX7.4C Multifunctional Platform Scale  
 HX7 HR Multifunctional Scale  
 XA 5Y.M.A.P Microbalance for pipettes calibration  
 WPY Multifunctional Scale  
 MYA 5Y.P Microbalance for pipettes calibration  
 C32 Multifunctional Scale  
 HY10.HRP.H High Resolution Scale  
 MA 5Y Moisture Analyzer  
 5Y.PM Precision balance  
 MYA 5Y.F Microbalance for Filters  
 HX7.4P2 H Multifunctional Stainless Steel Beam Scales  
 TMX H Stainless steel Multifunctional Scale  
 WLY Scales for PGC  
 AS X7 Analytical Balance  
 PS X7 Balances for PGC  
 TMX Multifunctional Scale  
 XA 5Y.M.A Microbalance  
 PS R2 Precision Balance  
 WLC X7 Precision balance  
 HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale  
 CY10 Scales for PGC  
 DWT/RC/HYF Checkweighers  
 PUE CY10 Weighing Terminal  
 CY10 Multifunctional Scale  
 HX7.4P H Multifunctional Stainless Steel Pallet Scales  
 HY10 Multifunctional Scale  
 HY10 H Stainless Steel Multifunctional Scale  
 HY10.HRP High Resolution Scale  
 C32.PM Precision Balance  
 PUE 7.1.HRP High Resolution Scale  
 Automatic Balances for. Pipettes Calibration  
 HX7 Multifunctional Scale  
 MA X7.A Moisture Analyzer  
 HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version  
 PS X7 Precision Balance  
 WPY Scales for PGC  
 XA 5Y.A Analytical Balance  
 HX7.4 H Multifunctional Stainless Steel Platform Scale

UYA 5Y Ultra-Microbalance  
 MA 5Y.IC Moisture Analyzer  
 MA R Moisture Analyzer  
 AS 5Y Analytical Balance  
 AK-4 Automatic Mass Comparator  
 PS 5Y Balances for PGC  
 AKM-2 Automatic Mass Comparator  
 PS R2.H Precision Balance  
 HY10.H scales for PGC  
 CY10 Precision balance  
 Analytical Balance for Filters XA 5Y.F  
 WLY Precision Balance  
 UMA Automatic Mass Comparator  
 APP 5Y.KO Manual Mass Comparator  
 Robotic Mass Comparators  
 Density measurement  
 MYA 5Y Microbalance  
 PS R1 Precision Balance  
 XA 5Y Analytical Balance  
 PUE HX7 Weighing Terminal  
 5Y.PM.KB Manual Mass Comparator  
 AS R2 PLUS Analytical Balance  
 PUE C32 Weighing Terminal  
 PUE 7.1 Weighing Terminal  
 PUE 5 Weighing Terminal  
 HRP 5Y.KO Mass Comparator  
 PS 5Y Precision Balance  
 Automatic Vacuum Mass Comparators  
 PMV 5Y Moisture Analyzer  
 XA 5Y.KO Manual Mass Comparator  
 HRP 5Y.KB Mass Comparator  
 XA 5Y.M Microbalance  
 UYA 5Y.KO Manual Mass Comparator  
 AS R1 PLUS Analytical Balance  
 UYA 5Y.F Ultra-Microbalance for Filters  
 H315 H Waterproof Platform Scale  
 HX7.4P Multifunctional pallet scales  
 Nano-Comparator NANO.AK-4.500.5Y  
 XA 5Y.M.A.S Microbalance for stents weighing