

LS2208 Barcode Scanner

WX-006-0005





The drawings, photos and graphics used are for illustrative purposes only.

Overview

LS2208 barcode scanner is a manual device designed for operation in adverse industrial conditions. It integrates with RADWAG terminals, PUE H315, PUE 5, PUE 7.1, HY10, PUE C32, and balances, 5Y, R2, X2, H7.

LS2208 is a high quality handheld barcode scanner. It has been designed for intense and continuous operation in retail sales, libraries, health centers and pharmacy, etc. Its light and ergonomic casing provides comfort of use. The function of dynamic scanning and wide scanning range, enable LS2208 to collect barcodes in a quick and reliable way, increasing work pace and efficiency.Intellistand™ system provides scanning of codes placed in front of the scanner, and easy switching between handheld and handfree mode (operation on a stand). A patented design of a single printed-circuit board eliminates failure endangered parts, i.e. connections between multiple boards. Scanning part of LS2208 with liquid polymere protects its module from shock based damages. LS2208 scanner features like integrated operation of various systems, universal cable or new symbols detection guarantee indispensable protection of an investment now and in the future.

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer MUYA Modules

HX7.4P2 Multifunctional Beam Scales WAY 5Y.KO Manual Mass Comparator

HX7 H Multifunctional Scale

HX7.4N H Multifunctional Stainless Steel Ramp Scale

HX7.4C Multifunctional Platform Scale

HX7 HR Multifunctional Scale

XA 5Y.M.A.P Microbalance for pipettes calibration

WPY Multifunctional Scale

MYA 5Y.P Microbalance for pipettes calibration

C32 Multifunctional Scale

HY10.HRP.H High Resolution Scale

MA 5Y Moisture Analyzer 5Y.PM Precision balance

MYA 5Y.F Microbalance for Filters

HX7.4P2 H Multifunctional Stainless Steel Beam Scales

TMX H Stainless steel Multifunctional Scale

WLY Scales for PGC
AS X7 Analytical Balance
PS X7 Balances for PGC
TMX Multifunctional Scale
XA 5Y.M.A Microbalance
PS R2 Precision Balance
WLC X7 Precision balance

HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale

CY10 Scales for PGC

DWT/RC/HYF Checkweighers PUE CY10 Weighing Terminal CY10 Multifunctional Scale

HX7.4P H Multifunctional Stainless Steel Pallet Scales

HY10 Multifunctional Scale

HY10 H Stainless Steel Multifunctional Scale

HY10.HRP High Resolution Scale C32.PM Precision Balance

PUE 7.1.HRP High Resolution Scale

Automatic Balances for. Pipettes Calibration

HX7 Multifunctional Scale MA X7.A Moisture Analyzer

HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version H315 H Waterproof Platform Scale

PS X7 Precision Balance WPY Scales for PGC XA 5Y.A Analytical Balance

HX7.4 H Multifunctional Stainless Steel Platform Scale

UYA 5Y Ultra-Microbalance MA 5Y.IC Moisture Analyzer MA R Moisture Analyzer AS 5Y Analytical Balance

AK-4 Automatic Mass Comparator

PS 5Y Balances for PGC

AKM-2 Automatic Mass Comparator

PS R2.H Precision Balance HY10.H scales for PGC CY10 Precision balance

Analytical Balance for Filters XA 5Y.F

WLY Precision Balance

UMA Automatic Mass Comparator APP 5Y.KO Manual Mass Comparator

Robotic Mass Comparators
Density measurement
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
PUE HX7 Weighing Terminal
5Y.PM.KB Manual Mass Comparator
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PUE 5 Weighing Terminal
HRP 5Y.KO Mass Comparator
PS 5Y Precision Balance

Automatic Vacuum Mass Comparators

PMV 5Y Moisture Analyzer

XA 5Y.KO Manual Mass Comparator HRP 5Y.KB Mass Comparator XA 5Y.M Microbalance

UYA 5Y.KO Manual Mass Comparator AS R1 PLUS Analytical Balance UYA 5Y.F Ultra-Microbalance for Filters H315 H Waterproof Platform Scale HX7.4P Multifunctional pallet scales

Nano-Comparator NANO.AK-4.500.5Y XA 5Y.M.A.S Microbalance for stents weighing