



More information on the website
mirror.radwag.com/de/info,w1,BUE

Prüfgewicht der E1 Klasse - Knopfgewicht - 50 g

OP-503-0051



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

| Messtechnische Parameter | |
|--------------------------|-----------|
| OIML-Klasse | E1 |
| Konstruktion | |
| Material | Edelstahl |
| Nominalwert | 50 g |

Kompatibel mit (Additional Fee)

Pinzetten

Pinself