



More information on the website  
[mirror.radwag.com/it/info,w1,BVY](https://mirror.radwag.com/it/info,w1,BVY)

# Cilindrico, set (1 mg - 500 g), cassa di legno

OK-501-0040



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Classe OIML	E1
Costruzione	
Materiale	Stainless steel
Peso	1 111,11 g
Valore nominale	1 mg - 500 g
Numero	24 pcs

## Compatibile con (Additional Fee)

Pinzette

Pennelli