



More information on the website
mirror.radwag.com/it/info,w1,C69

Cilindrico, set (1 mg - 2 kg), imballaggi in plastica

OK-102-0044



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri metrologici	
Classe OIML	E2
Costruzione	
Materiale	Stainless steel
Peso	6 111,11 g
Valore nominale	1 mg – 2 kg
Numero	27 pcs

Compatibile con (Additional Fee)

Pinzette

Pennelli