



More information on the website  
[mirror.radwag.com/it/info,w1,C7Q](https://mirror.radwag.com/it/info,w1,C7Q)

# Peso standard E1 - filo di acciaio inossidabile - 500 mg

OP-503-0036



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Classe OIML	E1
Costruzione	
Materiale	Stainless steel
Valore nominale	500 mg

## Compatibile con (Additional Fee)

Pinzette

Pennelli