



More information on the website
mirror.radwag.com/it/info,w1,C80

Pipetta automatica 10 μ l ÷ 100 μ l - RW8-104-20-9

PA-203-0004



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri metrologici	
Volume	10 – 100 μ l
Systematic error	3,0% / 0,8%
Random error	1,5% / 0,15%
Tips	300 μ l



Compatibile con (Additional Fee)

Microbilancia per la calibrazione delle pipette XA 5Y.M.A.P
 Microbilancia per la calibrazione delle pipette MYA 5Y.P

Bilance automatiche per la calibrazione delle pipette
 Puntali per pipette