



More information on the website  
[mirror.radwag.com/tr/info,w1,C80](https://mirror.radwag.com/tr/info,w1,C80)

# 10 µl ÷ 100 µl - RW8-104-20-9 Otomatik Pipetler

PA-203-0004



The drawings, photos and graphics used are for illustrative purposes only.

## Teknik Veriler

### Metrolojik parametreler

Hacim	10 – 100 µl
Systematic error	3,0% / 0,8%
Tekrarlanabilirlik error	1,5% / 0,15%
Uçlar	300 µl



## İle uyumlu (Additional Fee)

XA 5Y.M.A.P pipetlerini kalibre etmek için Mikro Terazi  
MYA 5Y.P pipetlerini kalibre etmek için Mikro Terazi

Pipet kalibrasyonu için otomatik teraziler  
Pipet uçları