



More information on the website  
[mirror.radwag.com/en/info,w1,CAI](http://mirror.radwag.com/en/info,w1,CAI)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 mg - 100 g), Wooden Box

OK-201-0072



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	211.11 g
Nominal value	1 mg – 100 g
Quantity	21 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes