

SDKP Workstation for Pipettes Calibration

WL-101-018²





The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters	
Delivery components	table for calibration, anti-vibration table, probe for temperature, humidity and pressure measurement, probe for water temperature measurement, sink for used water, armrests, LCD monitor, PC with pre-installed software PIPETTES RADWAG, power supply, keyboard and mouse, container for used water
Packaging dimensions W x D x H	820×1160×890 mm
Length	1100 mm
Width	750 mm
Height	1100 mm
Number of workstations	1

Compatible with (Additional Fee)

AS 5Y Analytical Balance AS R2 PLUS Analytical Balance AS R1 PLUS Analytical Balance THBR 2.0 System - Ambient Conditions Monitoring

Device dimensions W x D x H

