



More information on the website  
[mirror.radwag.com/en/info,w1,CBG](https://mirror.radwag.com/en/info,w1,CBG)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 100 g), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 211.11 g        |
| Nominal value           | 1 mg – 100 g    |
| Quantity                | 21 pcs          |

## Compatible with

Tweezers

Brushes