



More information on the website  
[mirror.radwag.com/en/info,w1,CBP](https://mirror.radwag.com/en/info,w1,CBP)

# F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 2 kg), Wooden Box

OK-201-0048



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 6 110 g         |
| Nominal value           | 1 g – 2 kg      |
| Quantity                | 15 pcs          |

## Compatible with (Additional Fee)

Tweezers

Brushes