



More information on the website  
[mirror.radwag.com/en/info,w1,D5Z](https://mirror.radwag.com/en/info,w1,D5Z)

# Knob weights, Set (1 mg - 200 g), Wooden Box

OK-501-0039



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class E1

### Construction

Material Stainless steel

Weight 611.11 g

Nominal value 1 mg – 200 g

Quantity 23 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes